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
FF – Filtering Facepiece
R – Resuable
NR – Non-Reusable
D – Dolomitt tested, clogging

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



Partikkelfiltre – mekanisme – Electret non-woven

Influence of Humidity on FFP Masks with Electret Filter Media Under Real-Life Wearing Conditions

Daniel Stoll and Sergiy Antonyuk


Figure 3. SEM images of the individual layers of a type A FFP2 mask.




Figure 1. Single fiber filtration mechanism.

Liu, C.; Dai, Z.; He, B.; Ke, Q.-F. The Effect of Temperature and Humidity on the Filtration Performance of Electret Melt-Blown Nonwovens. *Materials* **2020**, *13*, 4774. https://doi.org/10.3390/ma13214774

EKSEMPEL - N95 Surgical Respirator Mask Specifications


https://www.detmoldmedical.com/

EMW filtertechnik www.emw.de

Sieving Effects


Fanges opp

EMW filtertechnik www.emw.de

Inertial Effects


Treghet

EMW filtertechnik www.emw.de

Interception

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

EMW filtertechnik www.emw.de

Diffusion Effects


Diffusjon


EMW filtertechnik www.emw.de


EMW filtertechnik www.emw.de


EMW filtertechnik www.emw.de

EMW filtertechnik www.emw.de

Filtereffektivitet. Minimumskrav i standarden EN 143 (EN 149)

Kritiske partikkelstørrelsen:

P1: 80 % filtreringseffektivitet (80%)

P2: 94 % filtreringseffektivitet (94%)

P3: 99,95 % filtreringseffektivitet (99%)

Åndedrettsvern i andre land

- USA: NIOSH
 - Sju filterklasser
 - Tallet angir minimum filtreringsnivå
 - 95 minst 95%
 - 99 minst 99%
 - 100 minst 99,97%
 - N – not oil resistant
 - R – somewhat resistant to oil (8 timer)
 - P – strongly resistant to oil (40 timer eller 30 dager)
 - Surgical – fluid resistant to a certified level measured against a stream of artifical blood.

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Table 7 — Particle and combination filter classes

Filter class	Efficiency %	Penetration %
F5	≥99,99	≤0,01
F4	≥99,90	≤0,1
F3	≥99,00	≤1,0
F2	≥95,00	≤5,0
F1	≥80,00	≤20,0

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WIL

Research Article

The Effects of Temperature and Humidity on Electrostatic Changes in Respirators and Their Filtration Efficiency

Jimin Kim,[1] **Jongmin Park**,[1] **Perng-Jy Tsai**,[2] **and Chungsik Yoon**[1,3]

The filtering mechanism of respirators involves static electricity and physical mechanisms, such as inertial impaction, interception, and diffusion. Static electricity is used to reduce pressure resistance, but electrostatic effects and the environmental variables that affect them are not well studied. Therefore, this study evaluated the contribution of static electricity to the filtration efficiency of respirators and the effects of changes in temperature and humidity on static electricity and the resulting filtration efficiency. Eight respirators from four manufacturers—four for workers (1st class respirators) and four for citizens (KF-94 respirators)—were selected, all of which met legal standards in Korea. The filtration efficiency and pressure drop were evaluated at baseline and after removing static electricity with isopropanol. To study the effects of humidity and temperature on the efficiency of the respirators, the respirators were exposed to 30, 50, and 98% humidity while the temperature was fixed at 25°C and to temperatures of −30, 25, and 70°C while the humidity was fixed at 50%. Static electricity was measured using a surface-potential meter, and the filtration efficiency was measured using NaCl and paraffin oil. When static electricity was removed from the 1st class and KF-94 respirators, the filtration efficiency decreased by 21.7% and 19.5%, respectively. Humidity affected the change in static electricity slightly more than temperature, but since not all static electricity was removed, the filtration efficiency was not affected greatly. The changes in static electricity and filtration efficiency due to humidity were significant, but this was because the standard deviation was small. Even under 98% humidity, the filtration efficiency was only 1–2% lower than that at 30% humidity and still met the 94% filtration efficiency criterion for respirators.

atmosphere

MDPI

Article

Influence of Humidity on FFP Masks with Electret Filter Media Under Real-Life Wearing Conditions

Daniel Stoll * and Sergiy Antonyuk

The results of the study demonstrate that the filtration efficiency of the masks is susceptible to alterations in real-life wearing conditions and external influences, such as those mimicked by simulated coughing. In particular, the experimental results indicate that adherence to the recommended wearing times is crucial, as a brief drying period can enhance the protective efficiency of the masks. The length of this drying phase can vary greatly, depending on the material composition of the mask.

Type A Type B Type C

Figure 2. Different types of FFP2 masks studied.


Figure 6. Filtration efficiency for sealed FFP2 masks over time during simulated breathing at constant 98% RH, 30 °C, and 15 L·min^{-1}.


Figure 5. Filtration efficiency for a normal fit of the FFP2 masks over time during simulated breathing at constant 98% RH, 30 °C, and 15 L·min^{-1}.

materials

MDPI

Article

The Effect of Temperature and Humidity on the Filtration Performance of Electret Melt-Blown Nonwovens

Chao Liu [1,2,*], Zijian Dai [1], Bin He [2] and Qin-Fei Ke [1,*]

4. Conclusions

In summary, high temperature and high humidity treatment had no obvious effect on the fiber morphology and structure of the electret melt-blown material, so its pressure drop changed little. However, higher temperatures (90 and 110 °C) could accelerate the charge decay of the electret melt-blown material compared with the lower temperature, resulting in a significant decrease in its filtration efficiency. Besides, under the same temperature conditions, high humidity could accelerate the release of charge, so the filtration efficiency dropped more. After the sample was treated at 110 °C and 90% RH for 24 h, the filtration efficiency dropped from 95.49% to 38.16% at a flow rate of 14.16 cm s^{-1}, the QF is only 0.009 Pa^{-1}, and the surface potential decreased to the lowest value of −1.01 kV. Thus, during the process of production, transportation and storage of electret melt-blown nonwovens and its products, high temperature and high humidity conditions should be avoided to reduce charge decay and ensure that the final products can provide effective protection for users.




Figure 8. Surface potential decay of melt-blown samples treated by different conditions.

Filtration efficiency (%)
Treatment time (h)
50°C
70°C
90°C
110°C
50°C-90%
70°C-90%
90°C-90%
110°C-90%
Air flow rate : 5.33 cm s-1

Filtration efficiency (%)
Treatment time (h)
50°C
70°C
90°C
110°C
50°C-90%
70°C-90%
90°C-90%
110°C-90%
Air flow rate : 14.16 cm s-1




Figure 1. Single fiber filtration mechanism.

Liu, C.; Dai, Z.; He, B.; Ke, Q.-F. The Effect of Temperature and Humidity on the Filtration Performance of Electret Melt-Blown Nonwovens. *Materials* **2020**, *13*, 4774. https://doi.org/10.3390/ma13214774

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

Dräger
Filter Selection Guide
Dräger. Technology for Life®

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Kilde: 3M

Table 4: Examples of contaminants, their OELs (here: AGWs, valid in Germany) and filter recommendations

This is only a small choice of contaminants as example. For more information and a wider choice of contaminants please try our Dräger VOICE database of hazardous substances on the internet (www.draeger.com/voice).


Contaminants	OEL		Filter	Colour
	ppm	mg/m³	type	code
Acetaldehyde	50	91	AX (P3)	
Acetamide	–	–	A P3	
Acetic Acid	10	25	B [E] (P2)	
Acetic anhydride	5	21	A (P2)	
Acetone	500	1200	AX (P3)	
Acetoncyanohydrine	–	–	A (P2)	
Acetonitrile	20	34	A (P3)	
Acetyl chloride	–	–	B P2	
Acrolein	0.1	0.25	AX (P3)	
Acrylamide	–	–	A P3	
Acrylic acid	10	30	A (P2)	
Acrylnitrile	carcinogen	(cat. 2)	A (P3)	
Aldrin	–	0.25 E	A P3	
Allyl chloride	carcinogen	(cat. 3B)	AX (P3)	
1-Allyloxy-2, 3-epoxypropane	–	–	A (P2)	
Allylpropyldisulfide	2	12	B (P2)	

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







MSA Response Guide
webapps2.msasafety.com/ResponseGuide/Home.aspx
MSA
The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Location: Norway (English) (EN) . Click here to change your selection
Chemical Database, Respiratory Protection and Gas Detection Selection
Cartridge Life Expectancy Calculator
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |


MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,54 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:20 °C Humidity:40%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Step 5
Breakthrough Concentration
Results
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Contaminant
Name or CAS #
Concentration
ppm
mg/m3
TLV (Optional)
ppm
mg/m3
Add Chemical
Can't find a chemical? Contact Us.
Contaminant
Concentration
TLV
Gas Filter
Toluene
54 ppm
25 ppm
A
Remove
Note: The Maximum Use Concentration (MUC) is calculated by multiplying the Assigned Protection Factors (as given in EN 529) by the local Threshold Limit Value (MUC=APFxTLV).
Consequently, calculated MUC can be different according to the country selection.
Whenever those data are not available for the respective countries, the lowest value will prevail in the calculation.
«Back
Next »
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |


MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,261 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:20 °C Humidity:40%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Step 5
Breakthrough Concentration
Results
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Environmental Data
Temperature:
20
°C
°F
Humidity:
40
%
Atmospheric Pressure or Altitude:
Pressure Unit mm Hg hPa
Altitude Unit feet meters
760 mmHg (0 (sea level) - 1000 m)
585 mmHg (1001-1800 m)
492 mmHg (1801-3000 m)
Custom Value
«Back
Next »
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |

MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,54 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:20 °C Humidity:40%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Half Mask Standard EN 148-1 thread
90 A2B2
Step 5
Breakthrough Concentration
Average Breathing Rate 30 lpm
Breakthrough Concentration 10 % of TLV
Results
Estimated breakthrough time is:
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Select Mask and Cartridge
Choose a Mask Type:
Half mask
Choose a Cartridge:
90 A2B2
«Back
Next »
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |

MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,54 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:20 °C Humidity:40%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Half Mask Standard EN 148-1 thread 90 A2B2
Step 5
Breakthrough Concentration
Average Breathing Rate 60 lpm
Breakthrough Concentration 10 % of TLV
Results
Estimated breakthrough time is:
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Breakthrough Concentration
Average Breathing Rate
Low (30 lpm)
Moderate (60 lpm)
High (85 lpm and above)
Breakthrough Concentration
10% of TLV
50% of TLV
Custom
«Back
Calculate
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |

MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,54 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:20 °C Humidity:40%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Half Mask Standard EN 148-1 thread
93 A2B2 P3
Step 5
Breakthrough Concentration
Average Breathing Rate 60 lpm
Breakthrough Concentration 10 % of TLV
Results
Estimated breakthrough time is:
47 hours and 56 minutes
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Service Life Calculation Results
Toluene
Estimated breakthrough time is:
47 hours and 56 minutes
Email
Save
See the full Terms and Conditions
• MSA recommends replacing chemical cartridges every shift or sooner if the calculated service time exceed 8h (480 min)
• European Standard EN 529:2005 (Recommendations for selection, use, care and maintenance) gives guidance for selection and use of respiratory protective devices. In any case, applicable National Regulations must be followed, if existing.
• A filtering device should not be used in the following conditions:
- in confined spaces
- in case of oxygen deficiency
- when the locally required Maximum Use Concentration or IDLH value is exceeded
Back
« Start Over
The PDF pop-up has been blocked. Please check the pop-up settings at the top of your browser.
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |

MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,54 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:30 °C Humidity:40%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Half Mask Standard EN 148-1 thread
93 A2B2 P3
Step 5
Breakthrough Concentration
Average Breathing Rate 60 lpm
Breakthrough Concentration 10 % of TLV
Results
Estimated breakthrough time is:
46 hours and 42 minutes
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Service Life Calculation Results
Toluene
Estimated breakthrough time is:
46 hours and 42 minutes
Email
Save
See the full Terms and Conditions
• MSA recommends replacing chemical cartridges every shift or sooner if the calculated service time exceed 8h (480 min)
• European Standard EN 529:2005 (Recommendations for selection, use, care and maintenance) gives guidance for selection and use of respiratory protective devices. In any case, applicable National Regulations must be followed, if existing.
• A filtering device should not be used in the following conditions:
- in confined spaces
- in case of oxygen deficiency
- when the locally required Maximum Use Concentration or IDLH value is exceeded
Back
« Start Over
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |

MSA Response Guide
webapps2.msasafety.com/ResponseGuide/InternationalChemicalCalculator.aspx
MSA The Safety Company
Response® Guide
Chemical Database
Cartridge Life Expectancy Calculator
Contact Us
Step 1
Language and Regulation
Country: Norway
Standard: EN
Step 2
Contaminants, Concentrations, TLV
Toluene,54 ppm, 25 ppm
Step 3
Atmospheric Conditions
Temperature:20 °C Humidity:80%
Atmospheric Pressure or Altitude:
760 mm Hg
Step 4
Respirator and Cartridge Selection
Half Mask Standard EN 148-1 thread
93 A2B2 P3
Step 5
Breakthrough Concentration
Average Breathing Rate 60 lpm
Breakthrough Concentration 10 % of TLV
Results
Estimated breakthrough time is:
25 hours and 52 minutes
Cartridge Life Expectancy Calculator
Step 1
Step 2
Step 3
Step 4
Step 5
Results
Service Life Calculation Results
Toluene
Estimated breakthrough time is:
25 hours and 52 minutes
Email
Save
See the full Terms and Conditions
• MSA recommends replacing chemical cartridges every shift or sooner if the calculated service time exceed 8h (480 min)
• European Standard EN 529:2005 (Recommendations for selection, use, care and maintenance) gives guidance for selection and use of respiratory protective devices. In any case, applicable National Regulations must be followed, if existing.
• A filtering device should not be used in the following conditions:
- in confined spaces
- in case of oxygen deficiency
- when the locally required Maximum Use Concentration or IDLH value is exceeded
Back
« Start Over
© MSA 2016. All Rights Reserved.
Location: Norway (English) (EN) | MSA Website | Contact Us |

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EN 529:2005

3.6
nominal protection factor
number derived from the maximum percentage of total inward leakage permitted in relevant European Standards for a given class of respiratory protective device. The relationship between nominal protection factor and total inward leakage can be expressed as follows:

$$\text{nominal protection factor} = \frac{100}{\text{permitted maximum percentage total inward leakage}}$$

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Tabell - EN 529

Table C.1 — Nominal protection factors and assigned protection factors used in different countries

Standard	Description	Class	NPF	Assigned Protection Factors used in some countries				
				FIN	D	I	S	UK
EN 149	Filtering half mask	FF P1	4	4	4	4	4	4
		FF P2	12	10	10	10	10	10
		FF P3	50	20	30	30	20	20
[illegible]	[illegible]	[illegible]	[illegible]		[illegible]			[illegible]
EN 140 (Mask)	Half mask and quarter mask with filter	P1	4	4	4	4	4	4
		P2	12	10	10	10	10	10
Filters EN 141 *) EN 143 EN 371 *) EN 372 *) EN 14387 EN 12083		P3	48		30	30		20
		GasX	50	20	30	30	20	10
		GasX P1	4					
		GasX P2	12					
		GasX P3	48		30	-.-		10
[illegible]	[illegible]	[illegible]	[illegible]		[illegible]			[illegible]
EN 136 (Mask)	Full face mask (all classes)	P1	5	4	4	4	4	4
		P2	16	15	15	15	15	10
Filters EN 141 *) EN 143 EN 371 *) EN 372 *) EN 14387 EN 12083		P3	1000	500	400	400	500	40
		GasX	2 000	500	400	400	500	20
		GasX P1	5					
		GasX P2	16					
		GasX P3	1 000		400	-.-		20

Tabell - EN 529 forts.

Table C.1 — Nominal protection factors and assigned protection factors used in different countries

Standard	Description	Class	NPF	Assigned Protection Factors used in some countries				
				FIN	D	I	S	UK
EN 12941	Powered filtering device incorporating a hood or a helmet	TH1	10	5	5	5 [b]	5	10
		TH2	50	20	20	20 [b]	20	20
		TH3	500	200	100	200 [b]	200	40
EN 12942	Powered assisted filtering device incorporating full face mask, half mask or quarter mask	TM1	20	10	10	10 [b]	10	10
		TM2	200	100	100	100 [b]	100	20
		TM3	2 000	1 000	500	400 [b]	1 000	40
EN 14593-1	Compressed air line breathing		2 000	1 000	1 000	400	1 000	40

Nominell beskyttelesefaktor

- FF
 - P1 - 4
 - P2 - 12
 - P3 - 50
- TH
 - TH1P - 10
 - TH2P - 50
 - TH3P - 500
- TM
 - TM1P - 20
 - TM2P - 200
 - TM3P - 2000

- Halvmaske
 - P1 - 4
 - P2 - 12
 - P3 - 50
 - Gass - 50

- Helmaske
 - P1 - 5
 - P2 -16
 - P3 - 1000
 - Gass - 2000

Skjegg

Facial Hairstyles and Filtering Facepiece Respirators

Intended for workers who wear tight-fitting respirators

CLEAN SHAVEN ✓

STUBBLE ✗

LONG STUBBLE ✗

FULL BEARD ✗

FRENCH FORK ✗

DUCKTAIL ✗

VERDI ✗

GARIBALDI ✗

BANDHOLZ ✗

SOUL PATCH ✓

GOATEE ✓
(Careful! Chin hair may easily cross the seal)

CHIN CURTAIN ✗

EXTENDED GOATEE ✗

CIRCLE BEARD ✗

ANCHOR ✓ / BALBO ✓
(Careful! Chin hair may easily cross the seal)

VAN DYKE ✗

IMPERIAL ✗

SIDE WHISKERS ✓

MUTTON CHOPS ✗

HULIHEE ✗

HORSESHOE ✓
(Careful not to cross the seal)

ZAPPA ✓

WALRUS ✓

PAINTER'S BRUSH ✓

CHEVRON ✓

HANDLEBAR ✓

PENCIL ✓

LAMPSHADE ✓

ZORRO ✓

VILLAIN ✓
(Careful not to cross the seal)

WET NOODLE ✗

ENGLISH ✗

DALI ✗

Original image vector by fredrisher/Shutterstock.com

[1]If your respirator has an exhalation valve, some of these styles may interfere with the valve working properly if the facial hair comes in contact with it.
[1]This graphic may not include all types of facial hairstyles. For any style, hair should not cross under the respirator sealing surface.
Source: OSHA Respiratory Protection Standard
https://www.osha.gov/pls/oshaweb/owadisp.show_document?p_table=standards&p_id=12716
Further Reading: NIOSH Respirator Trusted-Source Webpage
https://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/respsource3fittest.html

Centers for Disease Control and Prevention
National Institute for Occupational Safety and Health


2017

HSE

Health and Safety Executive

The effect of wearer stubble on the protection given by Filtering Facepieces Class 3 (FFP3) and Half Masks

Prepared by the **Health and Safety Laboratory** for the Health and Safety Executive 2015

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Kulde

Annals of Work Exposures and Health, 2026, 70, wxag063
https://doi.org/10.1093/annweh/wxag063
Published: 2 September 2026
Original article

BOHS British Occupational Hygiene Society
The Chartered Society for Worker Health Protection
OXFORD

Facial cooling while using respiratory protective devices in the cold

Jenni Kaisto[1,2*], Kirsi Jussila[1,2], Arto Reiman[2], and Sirkka Rissanen[3]

[1]Comprehensive Security in Working Life, Finnish Institute of Occupational Health, Työterveyslaitos, Oulu 90032, Finland
[2]Industrial Engineering and Management, University of Oulu, P.O. Box 8000, Oulu 90014, Finland
[3]Occupational Safety, Finnish Institute of Occupational Health, Työterveyslaitos 90032, Finland
*Corresponding author: Email: jenni.kaisto@ttl.fi.


Abstract

The use of respiratory protective devices (RPDs) is essential in various work environments, including those with dust generation. When the work environment is cold, the low ambient temperature can accentuate discomfort from cooling and may influence workers' selection of RPDs, potentially affecting protective efficiency. This study aimed to quantify facial cooling while wearing particle-filtering facepiece respirators in cold conditions. Six different RPDs were selected: two filtering facepiece respirators (FFRs), two powered

What's Important About this Paper?

This study investigated facial cooling while using respiratory protective devices in cold conditions. The results showed that use of powered air-purifying respirators with continuous airflow at −20 °C causes excessive facial cooling and discomfort for the user. While selecting respiratory protective devices for use, the level of protection needed should always be prioritised but this study indicates that working temperature should also be taken into consideration.


Figure 1 Mean facial skin temperature while wearing different types of RPDs at −20 ° C during standing (0 to 15 and 30 ...

Ann Work Expo Health, Volume 70, Issue 7, October 2026, wxag063, https://doi.org/10.1093/annweh/wxag063

The content of this slide may be subject to copyright: please see the slide notes for details.


Table 2 Studied RPD types and their abbreviations and descriptions.

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
Kilde: 3M

Journal of Occupational and Environmental Hygiene, 8: 267–270
ISSN: 1545-9624 print / 1545-9632 online
Copyright © 2011 JOEH, LLC
DOI: 10.1080/15459624.2011.566016

Health Care Workers and Respiratory Protection: Is the User Seal Check a Surrogate for Respirator Fit-Testing?

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https://www.unn.no/ama/demaskus

https://doi.org/10.1016/j.shaw.2019.06.004

Contents lists available at ScienceDirect

Safety and Health at Work

journal homepage: www.e-shaw.org

Original Article

Do Various Respirator Models Fit the Workers in the Norwegian Smelting Industry?

Solveig Foereland [1,2,*], Oeystein Robertsen [3,4], Marit Noest Hegseth [4,5]

[1] Department of Occupational Medicine, St. Olavs Hospital, Trondheim University Hospital, Trondheim, Norway
[2] Department of Geoscience and Petroleum, NTNU, Norwegian University of Science and Technology, Trondheim, Norway
[3] Department of Psychology, UiT – the Artic University of Norway, Tromsø, Norway
[4] Department of Occupational and Environmental Medicine, University Hospital of North Norway, Tromsø, Norway
[5] Department of Community Medicine, UiT – the Artic University of Norway, Tromsø, Norway

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